

Protocol: Hemocytometer Cell Counting

Application:

Manual counting of cell suspensions


Procedure:

1. Aliquot 100 μL of 0.4% Trypan Blue solution into appropriately labeled tube.
2. Aliquot 100 μL of pre-mixed cell suspension into the appropriately labeled tube from Step 1.
 - a 1:1 ratio of Trypan Blue to cell suspension (equals 2-fold dilution)¹.
3. Vortex briefly (touch-spin if using 1.5 mL microfuge tubes to bring contents down from lid).
4. Pipette 10 μL of Trypan Blue-cell suspension into Hemocytometer chamber, between chamber notch and glass coverslip (**Figure 1A**, black arrow).
5. Count cells in boxes 1-4 (**Figure 1B**, red boxes)^{2,3}.
6. Calculate (A) cells/mL and (B) total cells in your sample using the formulas below:

$$(A) \text{ Cells/mL} = \left(\frac{\text{total cells counted}}{\text{number of boxes counted}} \right) \times \text{dilution factor} \times 10,000$$

$$(B) \text{ Total cells in sample} = \text{cells/mL} \times \text{total sample volume}$$

Example: $(100 \text{ cells counted} / 4 \text{ boxes}) \times 2 \times 10,000 = 500,000$ (or 5×10^5) cells/mL
 $5 \times 10^5 \text{ cells/mL} \times 10 \text{ mL sample volume} = 5,000,000$ (or 5×10^6) total cells in sample



¹ Trypan Blue-cell suspension dilution may change based upon predicted cell concentration.

² Aim to count ≥ 100 cells. If counting ≥ 100 cells/box, then less boxes can be counted assuming equal distribution of cells across the Hemocytometer chamber.

³ To avoid overestimation of cell numbers, include cells that lie on the top and right edges of each box. Do not include cells along the bottom and left edges.