

**Author:** Peter D Pioli  
**Last Edited:** 09/13/2019

## **Protocol: Transformation of DH5 $\alpha$ *Escherichia coli* Bacteria Cells**

### **Application:**

Transforming of recombinant DNA using DH5 $\alpha$  *E. coli* competent cells.

### **Procedure:**

1. Thaw competent cells on ice (source: Invitrogen MAX Efficiency DH5 $\alpha$ , cat# 18258012)
2. Aliquot 20  $\mu$ L of cells per 1.5 mL tube (negative control, positive control, experiments)
  - Refreeze unused cells on dry ice then into -80  $^{\circ}$ C.
3. Set up control tubes:
  - a) Negative control => no DNA
  - b) Positive control => 1  $\mu$ L of pUC19 DNA (0.01  $\mu$ g/ml stock)
4. Add 1  $\mu$ L of each new ligation (experiment tube) per DH5 $\alpha$  cell tube.
5. Finger flick tube bottom to mix
6. Ice for 30 min.
7. Heat shock 45 seconds at 42  $^{\circ}$ C in H<sub>2</sub>O bath.
8. Ice for 2 min.
9. Add 900  $\mu$ L **Luria-Bertani** (LB) media w/o Ampicillin.
10. Shake at 225 rpm at 37  $^{\circ}$ C for 1 hour.
11. For plating:
  - Positive and Negative controls: use 100  $\mu$ L of undiluted samples.
  - Each experimental tube/ligation, use:
    - a) 100  $\mu$ L undiluted (from above)
    - b) 100  $\mu$ L of 10x dilution (100  $\mu$ L undiluted + 900  $\mu$ L LB media)
12. Spread all conditions of DH5 $\alpha$  containing tubes onto LB + Amp (50  $\mu$ g/mL) plates.  
(negative control, positive control, 2x each experimental tube/ligation)  
**\*\*\*Performed near hot flame for sterility\*\*\***
13. Allow plates to dry on benchtop at room temperature.
14. Invert plates (agar side up) and incubate at 37  $^{\circ}$ C overnight.