

Protocol: Running Pooled 10x Genomics Libraries using Illumina MiSeq

- MiSeq Reagent Kit v3 supports 6-20 pM loading concentration. 75% of max capacity is recommended (~15 pM).
- Thaw the HT1 and Reagent Cartridge at 4 °C, o/n. Reagents are stable up to one week when stored at 4 °C if not used immediately.
- Start machine wash (3x) w/ MqH₂O + 0.5% Tween-20

Run ID: Pioli _____

	REF	LOT	Exp. Date
Box 1 of 2			
Box 2 of 2 (optional)			
FC ID:			
PhiX control			

- Dilution of NaOH: Always prepare *freshly diluted NaOH*. And to prevent small pipetting error, prepare at least 1 mL of freshly diluted NaOH (within 12 hrs).
 - ⇒ 1 N of NaOH stock: -20 °C Stock
 - ⇒ Need **0.2 N NaOH** (5x dilution) = 20 µL of 1 N NaOH + 80 µL of MqH₂O
 - ⇒ **Mix well.**
- Inspect the Reagent Cartridge:
 - Invert 10x to mix cartridge before use
 - Check the position 1, 2 and 4 to make sure they are fully mixed and no precipitates.
 - Tap the bottom of cartridge to get rid of bubbles if any.
- Re-boot the System.
- Making 4 nM Stock:

[Library 1] = _____ nM; need 10 µL of 4 nM => _____ µL + _____ µL of MqH₂O.

[Library 2] = _____ nM; need 10 µL of 4 nM => _____ µL + _____ µL of MqH₂O.

[PhiX] = 10 nM; need 5 µL of 4 nM => 2 µL + 3 µL MqH₂O.
- Denature a 4 nM **Pooled Libraries**:
 - 5 µL of 4 nM **Library 1** + 5 µL 4 nM **Library 2** + 10 µL of **0.2 N NaOH**
 - Vortex briefly and a quick spin (or ... centrifuge at 280 x g for 1 min).
 - Incubate at RT for 5 mins (exactly).

Note: Scale up accordingly if pooling >2 libraries. For example, pooling 3 libraries (5 μL of 4 nM each) would require 15 μL 0.2 N NaOH.

- Denature a 4 nM **PhiX**:
 - 5 μL of 4 nM **PhiX** + 5 μL of 0.2 N NaOH
 - Vortex briefly and a quick spin (or ... centrifuge at 280 x g for 1 min).
 - Incubate at RT for 5 mins (exactly).
- For a 20 pM **Library** or **PhiX**:
 - ⇒ Add 990 μL of prechilled HT1 to 10 μL of denatured **Pooled Libraries** or **PhiX**.
- Dilute Denatured **Library/PhiX** to 15 pM; 75% of max capacity

	10 pM	15 pM	18 pM	20 pM
20 pM library/PhiX	400 μL	600 μL	720 μL	800 μL
Prechilled HT1	400 μL	200 μL	80 μL	0 μL

Mix well, then quick spin.

- Combine **Library** and **PhiX**

	1% Spike-In	5% Spike-In	10% Spike-In	15% Spike-In
15 pM denatured/diluted PhiX (stored up to 3 wks at -20 °C)	7 μL	35 μL	70 μL	105 μL
15 pM denatured/diluted Library	693 μL	665 μL	630 μL	595 μL
Total	700 μL ; 600 μL needed. Leave on ice until loading to cartridge			

- Load sample library
 - Clean the position 17 “Load Sample” with a Kimwipe®.
 - Punch a hole on the foil seal with a clean 1-mL pipet tip.
 - Load 600 μL of **Sample-PhiX**. Avoid touching the foil.
- Setup “Manual” run
 - Run ID: Pioli_____
 - Sample Index read (i7): 8, NNNNNNNN
 - Paired-read
 - Read 1 (10x Barcode UMI): 26
 - Read 2 (Insert): 91
- Clean the Flow Cell
 - Remove the Flow Cell from its container.

- b. Rinse the Flow Cell with MqH₂O, **rinse well to get rid of excess salts!**
 - c. Be careful around the black Flow Cell port gasket!
 - d. Thoroughly dry with lint-free lens cleaning tissue.
 - e. (Optional) Clean the Flow Cell glass w/ alcohol wipe (not on the port gasket). Check the glass for streaks/fingerprints/lint/tissue fibers.
 - f. (Optional) Dry excess EtOH w/ lint-free lens cleaning tissue.
 - g. Visual check the cell port.
- Load the Flow Cell properly by holding the edges of Flow Cell. Wait for the RFID of Flow Cell identified! Close the Flow Cell compartment door.
 - Raise the sipper handle for loading the PR2 (at 4 °C) and check the Waste Bottle. Empty Waste Bottle if necessary, put back into position and then slowly lower the sipper handle. Wait for the RFID of the PR2 bottle.
 - Load the Reagent Cartridge (w/ sample) and close the chiller door. Wait for the RFID of the Reagent Cartridge.
 - Review the run parameters and perform a pre-run check before starting the run.
 - Start Run.
 - Perform a Post-Run Wash (1x) with MqH₂O + 0.5% Tween-20.